

Refeyn Mass Photometer standards

Standard	kDa	First dilution of stock	Dilution on 2MP gasket	Supplier
Beta-Amylase (sweet potato) 3 mg/ml	56, 112, 224	1 in 200	1 in 20 μ l	Sigma Cat no. A8781
Albumin, bovine serum 10 mg/ml	66	1 in 100	2 in 20 μ l	Merck Cat no. A8531
Urease	272, 544			SLS U7752
Thyroglobulin (bovine) 7 mg/ml	670, 1340	1 in 100	2 in 20 μ l	Sigma Cat no. T9145

Gaskets: 53 left

Slides: 60 left